Kentech Instruments Ltd.

Power supply for P860X Streak Tube

17th. March 2004

PLEASE READ THIS MANUAL CAREFULLY BEFORE USING THE SWEEPUNIT.

Serial Number J03****



Kentech Instruments Ltd.,

Unit 9, Hall Farm Workshops, South Moreton, Didcot, Oxon, OX11 9AG, U.K.

Tel: 01235 510 748 Fax: 01235 510 722

International tel: (44) 1235 510 748 International fax: (44)1235 510 722

E-mail akldb@kentech.co.uk

Contents

DISCI	DISCLAIMER	
1	INTRODUCTION	4
1.1	SPECIFICATIONS	4
2	CIRCUIT DESCRIPTIONS	4

Figure Captions

Figure 1	Main board components	5
Figure 2	Bottom tracks	6
Figure 3	Top Tracks	7
Figure 4	15kV Photocathode supply	8
Figure 5	B2 supply	9
Figure 6	F2 supply	10
Figure 7	QLV supply	11
Figure 8	QLH supply	12
Figure 9	Monitors on 5kV and 15kV rails	13
Figure 11	Main Power supply	14

DISCLAIMER

This equipment contains high voltage power supplies. Although the current supply capacity is small, careless use could result in electric shock. It is assumed that this highly specialised equipment will only be used by qualified personnel.

Kentech Instruments Ltd. accept no responsibility for any electric shock or injury arising from use or misuse of this equipment. It is the responsibility of the user to exercise care and common sense with this highly versatile equipment.

1 INTRODUCTION

This manual describes the operation and use of the Kentech Power supply for Photonis P860X Streak Tube. The manual gives the electrical specifications and describes the setting up procedure.

1.1 SPECIFICATIONS

Five output rails
Photocathode -15kV
Quadrapolar lens (QLV) +500V +/- 100V
Quadrapolar lens (QLH) -500V +/- 100V
Focus Electrode F2 -5kV +/- 500V
Focus Electrode B2 -5kV +/- 500V
Self latching HT enable
Two interlocks
universal mains power

2 CIRCUIT DESCRIPTIONS

The circuit has five separate high stability powers supplies. Each supply is controlled by a DC control voltage to set the votlage. In addition there are internal potentiometers that are used to set up the start and end voltage of the range available from the front panel, see figure??.

Zener diode chains are used to limit the output voltages in the event of a feedback loop failure of the stabalised supplies. In addition the supplies F2and B2 are linked with zener chains that limit the voltage between them to ??

Two interlocks are provided. Each output needs to be shorted to ground to allow the unit to turn on the HY. Interlock failure will turn of the HT. We do not recomend that these interlocks be relied upon to protect the image tube, but they will prevent accidental enabling of the supplies.

The supply rails, Photocathode, F2 and B2 also have currnt monitors on them. Excessive current drwn from these rails will cause the LED on the front panel to illuminateThe thresholds can be set inside the unit, see figure??

Each supply rail also has a built in low voltage monitor point hese have been calibrated to the voltage monitored at the output with a 1GW probe. Note that this may not be quite the voltage tat would appear on a tube as the tube will draw less current than the 1GW probe and there is a series current limiting resistor.


The monior outputs are calibrated on the assumption that they will be connected to a **WM**put impedance meter.

For rails Photocathode, F2 ad B2, the monitor output is 1/1000th of the output voltage. For rails QLV and QLH they are 1/100th.

The output of the photocathode rail can be adjusted via the front panel preset potentiometer with te use of a small screwdriverthe other fours rails have multiturn lockable knobs with counter They are setup so that the centre voltage (counter at 500) corresponds to the centre of the range as stated in the specification.


Kentech Instruments Ltd., Unit 9, Hall Farm Workshops, South Moreton, Didcot, Oxon, OX11 9AG, England. 17th. March 2004


4





Kentech Instruments Ltd., Unit 9, Hall Farm Workshops, South Moreton, Didcot, Oxon, OX11 9AG, England.


5 17th. March 2004

















12

